


Q1Scientific provides industry-leading temperature-controlled storage facilities and state-of-the-art technology solutions, to manage sample inventory. Q1 Scientific supports customers in maximising opportunities, minimising risk and reducing costs.



## Why use Q1Scientific

- Maximise space:** Space comes at a premium. Free up valuable space and resources in your factory to work on your core operations.
- Turn CAPEX to OPEX:** With a reduction in your capital investment for storage facilities and infrastructure.
- Reduce costs:** Our pay-per-use model eliminates waste, unseen maintenance, service and security costs with no unexpected facility costs.
- Unlimited capacity:** Scale up or down your stability storage with no minimum contract.
- Manage risk and compliance:** Through stability storage at our cGMP facilities with extensive validation, quality systems, 365/24/7 security monitoring, backup power and systems.
- Complete visibility:** Access your sample data 24/7 through the secure Q1 Scientific documentation portal.

## Other Services Provided

- Retain Sample Storage
- Clinical Trial Samples
- Biological Storage
- Accelerated and Stressed Conditions
- Freeze/Thaw Thermal Cycling
- Heritage & Conservation Storage
- Medical Device Storage
- Photostability Studies
- Disaster Recovery Services



*Your* **Trusted**  
*stability storage*  
**partner**

## New Stability Study

- Enquire by phone or email
- Unique conditions catered for

## Protocol

- Stability program to cater exactly to your needs
- Client has complete control
- Flexibility on length of study, conditions, number of pulls

## Sample Management Process

### Collection and Transportation to Q1 Scientific

- Q1 Scientific can collect samples from your site and transport to Q1.
- Samples can be placed on storage on same day.
- Client may arrange their own transport if required

### Samples on Stability

- Samples placed on stability as per client protocol
- Calendar of pull points logged as per clients protocol
- Complete sample visibility via client portal

### Delivery to Client Site/3rd Party Laboratory

- Delivery to client site or 3rd party laboratory for testing.
- Transport via Q1 Scientific trained operative or client collection
- Samples transported at frozen, chilled or ambient conditions

### Scheduled/ Unscheduled Sample Pull/Return

- Scheduled or unscheduled pulls per client calendar
- Confirmation requested from client before any sample removed

## Secure


- 24/7 monitored facility
- 100% network and power back-up
- Documented chain of custody

## Trusted

- Dedicated stability storage experts
- Validated equipment and processes
- GMP Quality Management System
- Web-based client portal for remote access and management

# Storage Conditions

Full ICH range of conditions. Customised conditions as required.



“Q1 Scientific offers state-of-the-art environmentally controlled and monitored storage facilities across the full range of temperatures and humidities. Stability storage conditions for all key ICH climatic zones; suitable for long term, intermediate and accelerated stability trials.”